

Alcon®

LADARVISION®4000

Excimer Laser System

FACTS YOU NEED TO KNOW ABOUT CustomCornea® LASER ASSISTED IN-SITU KERATOMILEUSIS (LASIK) SURGERY

PATIENT INFORMATION BOOKLET

For Nearsightedness (Myopia) up to -8.00D Sphere with -0.50D to -4.00D of
Astigmatism

Please read this entire booklet. Discuss its contents with your doctor so that you have all of your questions answered to your satisfaction. Ask any questions you may have before you agree to the surgery.

Alcon, Inc.
2501 Discovery Drive, Suite 500
Orlando, FL 32826 U.S.A.

Tel: (877) 523-2784
Fax: (407) 384-1677
www.ladarvision.com

Outside the U.S., contact your local Alcon office

Copyright 2004 by Alcon, Inc.

This booklet may be reproduced only by a treating physician for use with patients considering CustomCornea® Laser Assisted In-Situ Keratomileusis (LASIK) Surgery using the LADARVision®4000 system. All other rights reserved.

This page left blank intentionally.

TABLE OF CONTENTS

A. Introduction 5

B. How Does CustomCornea® LASIK Correct Nearsightedness (Myopia) with Astigmatism? 6

C. Benefits of CustomCornea® LASIK 8

 Clinical Study to Evaluate Benefits 8

 Patient Demographics 8

 Visual Acuity *without* Glasses After Surgery 8

 Visual Acuity *without* Glasses After Surgery and *with* Glasses Before Surgery . 9

D. Risks of CustomCornea® LASIK 9

 Contraindications – When Can’t You Have Surgery? 10

 What Warnings and Other Information Do You Need to Know About? 10

 Warnings 10

 Precautions 11

 During the First Week Following Surgery 12

 During One to Six Months Following Surgery 12

 Clinical Study to Evaluate Risks 12

 Visual Acuity *with* Glasses After Surgery 13

 Change in Visual Acuity *with* Glasses After Surgery 13

 Adverse Events and Complications 14

 Worse and Significantly Worse Symptoms After Surgery 15

E. Patient Questionnaire Responses 15

F. Are You A Good Candidate For CustomCornea® LASIK? 17

G. What Should You Expect During CustomCornea® LASIK Surgery? 17

 Before The Surgery 17

 The Day Of Surgery 17

 The First Days After Surgery 18

H. Questions To Ask Your Doctor 19

I. Self-Test 20

J. Summary Of Important Information 21

K. Glossary 22

L. Patient Assistance Information 25

M. Index 26

This page left blank intentionally.

A. Introduction

The purpose of this booklet is to provide you with information on laser eye surgery. Please read this entire booklet carefully. See the “Glossary” (Section K) for an explanation of words shown in *italics*. Discuss your questions with a doctor trained in laser eye surgery. You need to understand the benefits and risks of this surgery before making a decision to have surgery.

You may have nearsightedness if you have trouble seeing objects clearly when they are far away. *Nearsightedness*, which is also called *myopia*, is a type of condition that causes blurred vision. In addition to nearsightedness, you may have *astigmatism* if you see that parts of objects appear more blurred than other parts. Glasses, contact lenses, or eye surgery can correct nearsightedness with astigmatism to help you see distant objects more clearly.

The types of eye surgeries that are available to correct nearsightedness with astigmatism are *Radial Keratotomy (RK)*, *Photorefractive Keratectomy (PRK)*, and *Laser Assisted In-Situ Keratomileusis (LASIK)*. These surgeries may not meet the vision requirements for some careers, such as military service.

Eye surgery can help you see more clearly by changing the shape of the front surface of your *cornea*, which is the clear layer at the front of your eye. RK uses a scalpel to make fine cuts in the cornea. PRK and LASIK use a laser to reshape the cornea. For LASIK, an instrument called a *microkeratome* first cuts a thin flap of tissue from the front of your cornea. This *corneal flap* is folded back, and the laser removes tissue under the flap to change the shape of the front surface of your eye. Then the flap is put back in place for the eye to heal.

Your eyeglass prescription is the usual way to tell how much nearsightedness with astigmatism you have. Another way is to measure the shape of the *wavefront* of reflected light coming out of your eye. A wavefront measurement gives more information about your nearsightedness with astigmatism than an eyeglass prescription. A wavefront measures all of the *focusing errors* in your eye, including complex errors that eyeglasses cannot correct. These complex focusing errors are called “higher-order *aberrations*”.


Your doctor can use either your eyeglass prescription or a wavefront measurement to plan LASIK surgery. LASIK surgery based on the eyeglass prescription is called *Conventional LASIK*. LASIK surgery based on the wavefront is called wavefront-guided LASIK. *CustomCornea® LASIK* is wavefront-guided surgery with the LADARVision® 4000 Excimer Laser System.

LASIK surgery is permanent. You can have LASIK surgery on one eye at a time. The second eye may have surgery on the same day or later, depending upon your choice and your doctor’s advice. Discuss with your doctor whether you are a good candidate for CustomCornea® LASIK surgery.

B. How Does CustomCornea® LASIK Correct Nearsightedness (Myopia) with Astigmatism?

You see objects because your eye focuses light into images. Your eye works like a camera. The camera lens focuses light to form clear images onto film. Both the cornea and lens in the eye focus light onto the back surface of the eye, called the *retina*. Diagram 1 shows that distant vision is clear when light focuses correctly.

DIAGRAM 1: NORMAL EYE



**Light focuses on the retina.
Vision is clear.**

Nearsightedness is a type of focusing error that results in blurred distant vision. Light from a distant object focuses in front of the retina, rather than on the retina. Diagram 2 shows that distant vision is blurry when light focuses incorrectly in a nearsighted eye.

DIAGRAM 2: NEARSIGHTED EYE



Light focuses in front of the retina. Vision is blurred.

You may have nearsightedness combined with astigmatism, which is another type of focusing error that results in blurred distant and/or near vision. This condition occurs if the front of your eye is more curved in some directions than others. Light rays from an object focus at different points inside the eye so some parts of objects appear more blurred than other parts. For example, a person with astigmatism might confuse an “R” with a “P” or an “F” on a sign. Diagram 3 shows an example of how light rays focusing at different points in front of the retina may cause blurred vision in a nearsighted astigmatic eye.


DIAGRAM 3: ASTIGMATIC EYE



Light focuses at different points in front of the retina. Vision is blurred.

Wearing glasses and contact lenses help your eye focus light properly on the retina. LASIK surgery focuses light properly by reshaping the cornea. LASIK surgery uses an *excimer laser* to remove a tiny amount of tissue from the cornea. This type of laser does not change any other parts of the eye. Diagram 4 shows that distant vision is clearer after LASIK.

DIAGRAM 4: CORRECTION OF VISION AFTER LASIK



Light focuses on the retina after surgery. Vision is clear.

CustomCornea® LASIK uses a wavefront unique to your eye for treatment. This wavefront is used to guide the laser that reshapes the cornea to correct focusing errors. The doctor measures the wavefront by projecting light into your eye and measuring the reflected light that comes out of your eye.

The LADARVision® 4000 system uses a very small laser beam to reshape the cornea. To correct for nearsightedness with astigmatism, the cornea receives hundreds to thousands of laser pulses during LASIK surgery. The system must place the laser pulses accurately to precisely reshape the cornea. Precise shaping of the cornea requires tracking and compensating for eye movement during surgery.

Your eyes are constantly making small movements. Some of these movements are involuntary and you do not notice them. You cannot hold your eye perfectly still even if you try. The LADARVision® 4000 system tracks and adjusts for eye movement during surgery. A high-speed active eye tracking system, called the LADARTracker® system, measures the eye position 4000 times a second.

In a clinical study¹, eye movement during surgery using the LADARVision® 4000 system was evaluated for 554 eyes. This study showed that:

- All eyes moved during surgery.
- The LADARTracker® system adjusted for this eye movement. The results of the surgery were about the same for eyes with large or small eye movements.
- Active eye tracking with the LADARTracker® system improves the accuracy of corneal shaping.

Without a system to track eye movements, any movement of the eye during surgery could move it away from its correct position under the laser beam. Before each laser pulse, the LADARTracker® system calculates where the eye has moved since the last pulse and moves the beam in exactly the same way, so the laser beam hits the cornea in the same place as if the eye had not moved.

¹ LADARVision® System PRK Myopia and Astigmatism study

C. Benefits of CustomCornea® LASIK

CustomCornea® LASIK surgery can correct up to -8.00 *diopters* (D) of nearsightedness with -0.50D to -4.00D of astigmatism. If you have nearsightedness with astigmatism within this range, CustomCornea® LASIK surgery may help you clearly see distant objects without eyeglasses or contact lenses.

Clinical Study to Evaluate Benefits

A clinical study was done to evaluate the benefits and risks of CustomCornea® LASIK. The study included 232 eyes to determine benefits and 331 eyes to determine risks. The study results are shown below and in “Risks of CustomCornea® LASIK” (Section D).

Patient Demographics

Table 1 shows the age, race, gender, and contact lens history of patients in the study.

| Table 1. Demographics of 331 Eyes of 167 Study Patients | | | | | | | |
|---|----------------|-------------|-------|--------|-------|----------------------|-------|
| Age | | Race | | Gender | | Contact Lens History | |
| Average: | 37 ± 9 years | Asian | 2.4% | Female | 40.8% | None | 24.2% |
| Range: | 19 to 57 years | Black | 0.6% | Male | 59.2% | Hard | 10.6% |
| | | Caucasian | 94.6% | | | Soft | 65.3% |
| | | East Indian | 0.6% | | | | |
| | | Hispanic | 1.8% | | | | |

Visual Acuity *without* Glasses After Surgery

Visual acuity measures the sharpness of vision using a letter chart. Table 2 shows that at least 97% of patients treated for nearsightedness with astigmatism saw 20/40 or better **without** glasses after surgery. Most states require that your vision be 20/40 or better if you drive **without** any glasses or contact lenses.

| Table 2. Visual Acuity <i>without</i> Glasses After Surgery | | |
|---|---------------------|---------------------|
| % of Eyes With: | 3 Months (N=232) | 6 Months (N=232) |
| 20/20 or better* | 81% | 86% |
| 20/20 or better | 79% | 84% |
| 20/25 or better | 96% | 94% |
| 20/40 or better | 99% | 97% |

N is the number of eyes studied.

* if 20/20 or better with glasses or contact lenses before surgery (N= 225 eyes).

Visual Acuity *without* Glasses After Surgery and *with* Glasses Before Surgery

Table 3 shows that 60.3% of the patients at 3 months and 67.2% at 6 months saw as well or better **without** glasses after CustomCornea® surgery as **with** glasses before surgery.

| Table 3. Visual Acuity <i>without</i> Glasses After Surgery Compared to <i>with</i> Glasses Before Surgery | | |
|---|-----------------------------|-----------------------------|
| % of Eyes With: | 3 Months (N=232) | 6 Months (N=232) |
| 2 lines better vision than with glasses | 0.9% | 1.3% |
| 1 line better vision than with glasses | 15.9% | 21.6% |
| the same vision as with glasses | 43.5% | 44.4% |
| 1 line worse vision than with glasses | 26.3% | 23.7% |
| 2 lines worse vision than with glasses | 11.2% | 4.3% |
| more than 2 lines worse vision than with glasses | 2.2% | 4.7% |

N is the number of eyes studied.

D. Risks of CustomCornea® LASIK

If you are not satisfied with your surgery results, your doctor may suggest another surgery. No data are available for CustomCornea® LASIK retreatments.

IMPORTANT: You may need to wear glasses or contact lenses for some activities after surgery. CustomCornea® LASIK does not take away the need for reading glasses. You may need reading glasses after CustomCornea® LASIK even if you did not need them before.

In some cases, your best vision **with** your glasses or contact lenses may be worse after CustomCornea® LASIK surgery than it was before surgery.

A number of risks from LASIK surgery are related to the corneal flap rather than the laser treatment. Some specific problems include: cutting an incomplete or irregular flap, loss of the flap, misalignment of the flap, and cutting all the way through the cornea with the microkeratome. These problems can lead to other complications, such as infections, *cataracts*, and permanent scarring or deformity of the eye.

Contraindications – When Can't You Have Surgery?

If you have any of the following situations or conditions, the risk of LASIK surgery is greater than the benefit. You should **NOT** have LASIK surgery if you:

- are pregnant or nursing. These conditions may cause temporary and unpredictable changes in your cornea that may interfere with the accuracy of the measurement of your cornea before the LASIK procedure.
- have a *collagen vascular* (e.g., rheumatoid arthritis), *autoimmune* (e.g., lupus), or *immunodeficiency disease* (e.g., AIDS). These conditions affect your body's ability to heal and may result in inflammation or swelling of parts of the body such as muscles, joints, and blood vessels.
- show signs of *keratoconus* or any other condition that causes a thinning of your cornea. This unstable condition of the cornea makes it unsafe to do LASIK procedures on eyes with this condition.
- are taking medications with ocular side effects, such as isotretinoin (Accutane²) for acne treatment or amiodarone hydrochloride (Cordarone³) for normalizing heart rhythm. These medications may affect the accuracy of the LASIK procedure or the way your cornea heals after surgery. This may result in poor vision after surgery.

What Warnings and Other Information Do You Need to Know About?

Warnings

If you have any of the following conditions, you may have LASIK if your doctor evaluates the seriousness of your condition and believes the benefit of having LASIK is greater than the risk. Discuss with your doctor if you have:

- diabetes. Diabetes may interfere with the healing of the cornea after LASIK.
- a history of *herpes simplex* or *herpes zoster* infection that has affected your eyes. LASIK may be more risky for patients who have an active or previous herpes infection that has affected their eyes.
- significant dry eye that is unresponsive to treatment. LASIK may increase the dry eye condition, which may or may not go away. This dryness may delay healing of the flap or interfere with the surface of the eye after surgery.
- severe allergies. The medications taken for severe allergies may interfere with the ability of the eye to heal after LASIK.

You will need eye drops to enlarge your pupil to at least 7mm to 11mm before surgery so the tracking system can more easily follow your eye during surgery. This effect of eye drops is only temporary.

² *Accutane* is a registered trademark of Hoffmann-La Roche Inc.

³ *Cordarone* is a registered trademark of Sanofi-Synthelabo

Precautions

If you have any of the following conditions, you should discuss this with your doctor. The safety and effectiveness for LASIK have **NOT** been established in patients:

- with unstable or worsening nearsightedness with astigmatism. Eyes with unstable nearsightedness with astigmatism are unable to be correctly measured to determine the right amount of the vision correction to provide.
- with conditions that may interfere with the ability to properly measure the eye to determine the right amount of vision correction, and may also affect the healing of the eye after the surgery, such as:
 - disease or corneal condition (for example, scar, infection, etc.).
 - injury to the center of the cornea where LASIK will reshape the cornea.
 - previous surgery on the cornea or inside the eye (for example, cataract surgery).
 - prior history of surgery to correct vision (for example, RK, PRK, LASIK).
- with a cornea that is too thin for LASIK to be completed safely. A flap needs to be cut into the cornea for the LASIK procedure. A proper flap cannot be created on a thin cornea.
- with a history of *glaucoma* (a condition usually associated with high eye pressure with damage to the nerve in the eye and possible loss of vision). It is unknown whether LASIK is safe for eyes with glaucoma.
- who are taking the medication sumatriptan (Imitrex⁴) for migraine headaches. It is unknown whether the use of this medication will interfere with the accuracy of the measurement of your cornea prior to LASIK or the healing of the eye after LASIK.
- under 21 years of age because it is unknown if the eye has reached its adult vision refraction. This may result in measurement of the amount of correction to provide being incorrect.
- over the long term (more than 6 months).
- with greater than -8.00D sphere of nearsightedness with less than -0.50D of astigmatism or greater than -4.00D of astigmatism. Corrections falling outside of the approved range have not been studied.
- for retreatment with this laser for LASIK. Retreatments have not been done enough times to allow an understanding of whether it is safe and effective.

Let your doctor know if you are taking any prescription medicines or any medications you bought without a prescription. These medications may interfere with the measurement prior to LASIK or the healing of the eye after LASIK.

The safety and effectiveness of wavefront-guided LASIK have only been established with an optical zone of 6.5mm and a total treatment zone of 9mm.

Before surgery, your doctor should evaluate your pupil size under dim lighting conditions. If your pupils in dim light are greater than the optical zone (>6.5mm) proposed by your doctor, consult with your doctor about the risk that the surgery may cause negative effects on your vision, such as glare, halos, and night driving difficulty.

⁴ Imitrex Reg TM of Glaxo Group Limited

Your doctor should also evaluate you for dry eyes before surgery. You may have dry eyes after LASIK surgery even if you did not have dry eyes before surgery.

In the clinical study, contrast sensitivity was measured in daylight and in dim light to determine how well patients can see in poor contrast conditions such as very dim light, rain, snow, and fog. The results are shown below:

Under daylight testing conditions, contrast sensitivity after CustomCornea® LASIK surgery as compared to before surgery showed a:

- gain in 6.0% of patients at 3 months and 6.5% at 6 months.
- loss in 6.5% of patients at 3 months and 1.7% at 6 months.

Under dim light testing conditions, contrast sensitivity after CustomCornea® LASIK surgery as compared to before surgery showed a:

- gain in 18.5% of patients at 3 months and 22.4% at 6 months.
- loss in 8.6% of patients at 3 months and 7.3% at 6 months.

During the First Week Following Surgery

- You may feel pain, discomfort, or have a feeling that something is in your eye. It may last up to 7 days after surgery.
- Your vision may be blurry or you may become more sensitive to light as your eye heals.
- You may have temporary swelling of the front surface of your eye.
- The pressure inside your eye may increase, usually due to the use of *anti-inflammatory medication* (eye drops) after surgery. Using another medication or stopping the anti-inflammatory medication can control the abnormal increase in eye pressure.

During One to Six Months Following Surgery

- Your vision should be stable 3 months after surgery. Some patients may notice that their vision improves or worsens. These small changes may occur up to 3 months or more after surgery. You should contact your doctor if you notice any change or loss of vision.
- You may become more sensitive to light. You may notice glare or have difficulty in driving at night.
- You may experience some dryness.

Clinical Study to Evaluate Risks

In the clinical study on CustomCornea® LASIK, vision **without** glasses improved for all eyes. Some people still needed glasses or contact lenses after surgery. At 3 months after surgery, 7% of patients had greater than 1D and 0.4% of patients had greater than 2D of correction over the intended amount. Under testing conditions using eyedrops to reduce focusing ability, 18% of patients had greater than 1D and 2% of patients had greater than 2D of correction over the intended amount at 3 months.

Visual Acuity *with* Glasses After Surgery

Table 4 shows that all patients from the study saw 20/32 or better **with** glasses at 3 and 6 months after surgery. Of those patients who saw 20/20 or better **with** glasses or contact lenses before surgery, 98.8% of patients at 3 months and 99.4% of patients at 6 months saw 20/20 or better **with** glasses or contact lenses after surgery.

| % of Eyes With: | 3 Months (N=331) | 6 Months (N=331) |
|------------------------|-----------------------------|-----------------------------|
| 20/20 or better* | 98.8% | 99.4% |
| 20/20 or better | 97.6% | 98.8% |
| 20/25 or better | 100% | 99.7% |
| 20/32 or better | 100% | 100% |

N is the number of eyes studied.

* if 20/20 or better with glasses or contact lenses before surgery (N=323 eyes).

Change in Visual Acuity *with* Glasses After Surgery

Under dim room lighting conditions, the best vision **with** glasses was measured using a standard (high-contrast) visual acuity chart and a 10% *low contrast visual acuity* chart. A standard chart has black letters on a white background. A 10% low contrast visual acuity chart has gray letters on a white background. Black letters are easier to see than gray letters. Low contrast acuity testing is another way to determine how well patients can see in poor contrast conditions such as very dim light, rain, snow, and fog. Table 5 compares the change in vision **with** glasses at 3 and 6 months before and after surgery.

| % of Eyes With: | Standard Chart | | 10% Low Contrast Chart | |
|---------------------------|------------------------------|------------------------------|-------------------------------|------------------------------|
| | 3 Months (N= 331) | 6 Months (N= 331) | 3 Months (N= 331) | 6 Months (N= 331) |
| loss of more than 2 lines | 0.0% | 0.0% | 0.9% | 0.6% |
| loss of 2 lines | 0.0% | 0.0% | 4.2% | 3.0% |
| loss of 1 line | 7.3% | 8.5% | 14.8% | 8.8% |
| no change | 55.9% | 47.7% | 37.2% | 38.7% |
| gain of 1 line | 35.0% | 41.4% | 36.6% | 38.1% |
| gain of 2 lines | 1.8% | 2.4% | 6.0% | 8.8% |
| gain of more than 2 lines | 0.0% | 0.0% | 0.3% | 2.1% |

N is the number of eyes studied.

Adverse Events and Complications

Some patients from the clinical study experienced adverse events and complications after CustomCornea® LASIK surgery as shown in Table 6.

| Table 6. Adverse Events and Complications | |
|--|------|
| Greater than or equal to 1% of eyes (N=331) had: | |
| <i>Corneal swelling</i> between one week and less than one month | 1.5% |
| Less than 1% of eyes (N=331) had: | |
| Cells growing under the flap | 0.9% |
| <i>Inflammation</i> of the cornea under the flap | 0.9% |
| <i>Corneal swelling</i> of the flap at one month or later | 0.6% |
| Double or ghost images | 0.6% |
| Feeling of something in the eye at one month or later | 0.6% |
| <i>Corneal wrinkle</i> in the flap | 0.3% |
| Flap creation without a hinge with the microkeratome | 0.3% |
| Irregular flap creation with the microkeratome | 0.3% |

N is the number of eyes studied.

There were no reports of the following adverse events and complications in the clinical study:

- eye pain at one month or later;
- corneal scratch at one month or later;
- corneal infection;
- corneal cloudiness at six months with a loss of 2 or more lines of visual acuity with glasses;
- loss of more than 10 letters (more than 2 lines) of visual acuity with glasses at six months;
- cells growing under the flap with a loss of 2 or more lines of visual acuity with glasses;
- breakdown of the flap;
- misaligned flap;
- eye pressure more than 25 mmHg;
- increase in eye pressure of more than 10 mmHg compared to before surgery;
- separation of the retina from the back of the eye;
- blockage of blood vessels in the retina.

Worse and Significantly Worse Symptoms After Surgery

Patients who were treated for nearsightedness with astigmatism rated the change in the following symptoms after surgery **without** glasses or contact lenses as worse or significantly worse compared to before surgery **with** glasses or contact lenses (Table 7).

| Table 7. Symptoms <i>without</i> Glasses After Surgery Compared to <i>with</i> Glasses Before Surgery | | | | |
|---|------------------|---------------------|------------------|---------------------|
| | 3 Months (N=232) | | 6 Months (N=232) | |
| Symptom | Worse | Significantly Worse | Worse | Significantly Worse |
| Blurring of Vision | 12.1% | 1.3% | 15.9% | 0.4% |
| Burning | 5.2%* | 0.0%* | 3.5%** | 0.9%** |
| Double Vision | 5.2%** | 2.6%** | 6.9% | 0.4% |
| Dryness | 23.3% | 2.2% | 17.7% | 1.7% |
| Excessive Tearing | 0.0%* | 0.0%* | 0.0%* | 0.0%* |
| Fluctuation of Vision | 22.4% | 2.2% | 16.4% | 1.7% |
| Glare | 20.8%** | 1.7%** | 13.4% | 0.0% |
| Gritty Feeling | 9.5% | 0.9% | 4.3% | 0.0% |
| Halos § | 19.8% | 1.3% | 16.4% | 0.9% |
| Headache | 5.2% | 0.0% | 2.6% | 0.9% |
| Light Sensitivity | 21.6% | 0.9% | 15.9% | 0.4% |
| Night Driving Difficulty | 9.5% | 1.7% | 10.3% | 0.0% |
| Pain | 3.4% | 0.0% | 2.6%* | 0.0%* |
| Redness | 5.6% | 0.4% | 1.3% | 0.0% |

N is the number of eyes studied.

For the symptoms noted with asterisks, *230 patients responded and **231 responded.

§ *Halos* are circular flares or rings of light that may appear around a headlight or other lighted object.

E. Patient Questionnaire Responses

After treatment for nearsightedness with astigmatism, 94% of patients reported at 3 months and 95% reported at 6 months that they never wore glasses or contact lenses. Quality of vision **without** glasses or contact lenses after surgery was compared to **with** glasses or contact lenses before surgery. Quality of vision was rated as unchanged, better, or significantly better than before surgery in 93% of patients at 3 months and 94% at 6 months. Approximately 88% of patients at 3 and 6 months reported they were satisfied or extremely satisfied with their results.

Patients rated the change in the following symptoms after surgery **without** glasses or contact lenses compared to before surgery **with** glasses or contact lenses. Table 8 shows that more than half of all patients being treated for nearsightedness with astigmatism reported that their symptoms were the same at 3 and 6 months after CustomCornea® LASIK surgery **without** glasses as before surgery **with** glasses.

| Table 8. Symptoms <i>without</i> Glasses After Surgery Compared to <i>with</i> Glasses Before Surgery | | | | | |
|--|-----------------------------|---------------|------------------|--------------|----------------------------|
| Symptom | Significantly Better | Better | No Change | Worse | Significantly Worse |
| 3 Months (N=232) | | | | | |
| Blurring of Vision | 9.1% | 10.3% | 67.2% | 12.1% | 1.3% |
| Burning* | 6.1% | 11.7% | 77.0% | 5.2% | 0.0% |
| Double Vision** | 7.4% | 3.0% | 81.8% | 5.2% | 2.6% |
| Dryness | 7.8% | 14.7% | 52.2% | 23.3% | 2.2% |
| Excessive Tearing* | 5.2% | 6.1% | 88.7% | 0.0% | 0.0% |
| Fluctuation of Vision | 9.1% | 10.8% | 55.6% | 22.4% | 2.2% |
| Glare** | 6.1% | 14.3% | 57.1% | 20.8% | 1.7% |
| Gritty Feeling | 7.8% | 9.9% | 72.0% | 9.5% | 0.9% |
| Halos § | 6.5% | 14.2% | 58.2% | 19.8% | 1.3% |
| Headache | 6.0% | 9.9% | 78.9% | 5.2% | 0.0% |
| Light Sensitivity | 4.3% | 14.7% | 58.6% | 21.6% | 0.9% |
| Night Driving Difficulty | 7.3% | 21.6% | 59.9% | 9.5% | 1.7% |
| Pain | 6.0% | 8.6% | 81.9% | 3.4% | 0.0% |
| Redness | 6.5% | 10.3% | 77.2% | 5.6% | 0.4% |
| 6 Months (N=232) | | | | | |
| Blurring of Vision | 10.8% | 10.3% | 62.5% | 15.9% | 0.4% |
| Burning** | 7.4% | 11.3% | 77.1% | 3.5% | 0.9% |
| Double Vision | 9.5% | 4.3% | 78.9% | 6.9% | 0.4% |
| Dryness | 11.2% | 14.2% | 55.2% | 17.7% | 1.7% |
| Excessive Tearing* | 6.5% | 4.3% | 89.1% | 0.0% | 0.0% |
| Fluctuation of Vision | 11.6% | 9.9% | 60.3% | 16.4% | 1.7% |
| Glare | 9.9% | 18.5% | 58.2% | 13.4% | 0.0% |
| Gritty Feeling | 12.1% | 10.3% | 73.3% | 4.3% | 0.0% |
| Halos § | 10.8% | 17.7% | 54.3% | 16.4% | 0.9% |
| Headache | 9.9% | 12.1% | 74.6% | 2.6% | 0.9% |
| Light Sensitivity | 7.8% | 18.5% | 57.3% | 15.9% | 0.4% |
| Night Driving Difficulty | 15.9% | 20.3% | 53.4% | 10.3% | 0.0% |
| Pain* | 8.3% | 9.6% | 79.6% | 2.6% | 0.0% |
| Redness | 9.1% | 16.8% | 72.8% | 1.3% | 0.0% |

N is the number of eyes studied.

For the symptoms noted with asterisks, *230 patients responded and **231 responded.

§ *Halos* are circular flares or rings of light that may appear around a headlight or other lighted object.

F. Are You A Good Candidate For CustomCornea® LASIK?

If you are considering CustomCornea® LASIK, you must:

- be at least 21 years of age.
- have a healthy eye with no eye disease or corneal condition (for example, scar, infection, etc.).
- have up to -8.00D sphere of nearsightedness combined with -0.50D to -4.00D of astigmatism.
- have stable nearsightedness with astigmatism as documented by less than or equal to 0.50D change each year for at least one year before your eye examination before surgery.
- be able to lie flat on your back.
- be able to look at a blinking fixation light during the entire surgery.
- be able to have eye drops that numb your eye and enlarge your pupil.
- understand the risks and benefits of CustomCornea® LASIK compared to other available treatments for nearsightedness with astigmatism.
- be willing to sign an Informed Consent Form, if provided by your doctor.

G. What Should You Expect During CustomCornea® LASIK Surgery?

Before The Surgery

Before surgery, your doctor needs to determine your complete medical and eye history and check the health of both your eyes. As part of this exam, your doctor will use a computer program to map the front surface of your eye. This exam will determine if your eyes are healthy and if you are a good candidate for CustomCornea® LASIK.

WARNING: You must stop wearing any contact lenses at least 3 weeks before this eye examination. Failure to do this may affect surgical results.

Tell your doctor if you take any prescription and non-prescription medications or have any allergies. Ask your doctor if you should eat or drink right before the surgery. **You should also arrange for transportation since you must not drive right after the surgery.** Your doctor will let you know when your vision is good enough to drive again.

The Day Of Surgery

To prepare for surgery, your doctor will use the wavefront system to take a picture of your eye. This helps to determine where the laser should treat your cornea. Your doctor will put eye drops to dilate (enlarge) the pupil in your eye(s). After 30-40 minutes, your doctor will measure the wavefront unique to your eye to determine the amount of laser treatment you need.

Your doctor will then place numbing eye drops in the eye to be treated. Numbing drops are used to control pain during surgery. The effects of the numbing eye drops will wear off after about 45-60 minutes. Your doctor will ask you to lie on your back on the laser bed. The laser bed is a flat cushioned surface that can be moved to position you for surgery. Your doctor will instruct you to watch a blinking fixation light. Your doctor will place an instrument between your eyelids to hold them open during the surgery. A temporary shield will cover the eye that is not having surgery.

An instrument called the microkeratome creates a flap of tissue in the cornea. Then, your doctor will reposition your head and activate the LADARTracker® system to track your eye movement. Your doctor will ask you to look directly at the blinking light. The laser in the LADARVision®4000 system will remove small amounts of tissue from your cornea. During the laser treatment, you will hear a “clicking” sound of laser pulses. The tracking system will follow eye movements and allow the laser to continue the treatment. You will be under the laser for several minutes. The use of the laser will take about one minute. Overall, the surgery takes about 10 minutes.

IMPORTANT: You must continue looking at the blinking light throughout the treatment, even if your vision begins to become cloudy during the procedure.

After the surgery is complete, your doctor will place some eye drops in your eye. Your doctor may cover your eye with a *bandage contact lens* to help heal the eye. For your eye protection and comfort, your doctor may apply a patch or shield over your eye.

The First Days After Surgery

You may be mildly sensitive to light and have a feeling that something is in your eye. Sunglasses may make you more comfortable. Also, you may experience pain. Your doctor can prescribe pain medication to make you more comfortable during the first few days after the surgery. A plastic shield may be used to protect your eye after LASIK. You will need to use lubricants, *antibiotic*, and *anti-inflammatory medications* in the first few days.

IMPORTANT: Use the lubricants and eye medications as directed by your doctor. Your results depend upon you following your doctor’s instructions.

WARNING: Your doctor will monitor you for any side effects if you need to use a topical *steroid medication*. Possible side effects of prolonged topical steroid use are:

- *ocular hypertension* (an increase in the eye pressure);
- *glaucoma* (a condition usually associated with high eye pressure that results in damage to the nerve in the eye and possible loss of vision);
- *cataract formation* (an opacity or clouding of the lens inside the eye that can cause a loss of vision).

DO NOT rub your eyes for the first 3 to 5 days. Rubbing your eye may move the flap. If you notice any sudden decrease in your vision, you should contact your doctor immediately. The flap may have moved and the doctor may need to reposition the flap.

H. Questions To Ask Your Doctor

You may want to ask the following questions to help you decide if CustomCornea® LASIK with the LADARVision®4000 system is right for you:

- What are my other options to correct my nearsightedness with astigmatism?
- Will I have to limit my activities after surgery and for how long?
- What are the benefits of CustomCornea® LASIK for my amount of nearsightedness with astigmatism?
- What vision can I expect in the first few months after surgery?
- If CustomCornea® LASIK does not correct my vision, what is the possibility that my glasses would need to be stronger than before? Could my need for glasses increase over time?
- Will I be able to wear contact lenses after LASIK if I need them?
- How is LASIK likely to affect my need to wear glasses or contact lenses as I get older?
- Will my cornea heal differently if injured after having LASIK?
- Should I have LASIK surgery in my other eye?
- How long will I have to wait before I can have surgery on my other eye?
- What vision problems might I experience if I have LASIK only on one eye?
- Do I have significant dry eye or large pupils that could produce undesirable side effects after LASIK surgery?

Discuss the cost of surgery and follow-up care needs with your doctor. Most health insurance policies do not cover laser vision correction.

I. Self-Test

Are You An Informed And Educated Patient?

Take the test below to see if you can correctly answer the following questions after reading this booklet.

| | TRUE | FALSE |
|--|--------------------------|--------------------------|
| 1. LASIK surgery is risk-free. | <input type="checkbox"/> | <input type="checkbox"/> |
| 2. It does not matter if I wear my contact lenses before surgery when my doctor told me not to wear them. | <input type="checkbox"/> | <input type="checkbox"/> |
| 3. Since the LADARVision®4000 system tracks my eye movements, I do not have to fixate on the blinking light. | <input type="checkbox"/> | <input type="checkbox"/> |
| 4. After the surgery, there is a good chance that I will be less dependent on eyeglasses or contact lenses. | <input type="checkbox"/> | <input type="checkbox"/> |
| 5. I may need reading glasses after LASIK surgery, even if I did not need them before. | <input type="checkbox"/> | <input type="checkbox"/> |
| 6. There is a risk that I may lose some vision after LASIK surgery. | <input type="checkbox"/> | <input type="checkbox"/> |
| 7. It does not matter if I am pregnant. | <input type="checkbox"/> | <input type="checkbox"/> |
| 8. If I have an autoimmune disease, I am still a good candidate for LASIK surgery. | <input type="checkbox"/> | <input type="checkbox"/> |
| 9. Significant dry eye or large pupils may produce undesirable side effects after LASIK surgery. | <input type="checkbox"/> | <input type="checkbox"/> |

You can find the answers to Self-Test at the bottom of the next page.

J. Summary Of Important Information

- CustomCornea® LASIK is a permanent irreversible surgery to the cornea.
- You may need to wear glasses or contact lenses for some activities after surgery. CustomCornea® LASIK does not take away the need for reading glasses, even if you have never worn them before.
- Your vision must be stable before CustomCornea® LASIK surgery. You must provide written evidence that your nearsightedness with astigmatism has changed less than or equal to 0.50D each year for at least 1 year.
- Pregnant and nursing women should wait until they are not pregnant and not nursing to have CustomCornea® LASIK surgery.
- You would not be a good candidate if you have *autoimmune or collagen vascular diseases*. If you have a condition that makes wound healing difficult, you would not be a good candidate.
- CustomCornea® LASIK surgery has some risks. Please read and understand this entire booklet, especially the sections on Benefits and Risks before you agree to the surgery.
- Some other options to correct nearsightedness with astigmatism include glasses, contact lenses, RK, PRK, and Conventional LASIK.
- RK, PRK, Conventional LASIK or CustomCornea® LASIK may not meet the vision requirements of some occupations, such as military service.
- Before considering CustomCornea® LASIK surgery you should:
 - a. have a complete eye examination.
 - b. talk with at least one eye care professional about CustomCornea® LASIK, especially the potential benefits, risks, and complications. You should discuss the time needed for healing after CustomCornea® LASIK.

Answers to Self-Test Questions:

1. False (see Section D: Risks)
2. False (see Section G: Before the Surgery)
3. False (see Section G: The Day of Surgery)
4. True (see Section C: Benefits)
5. True (see Section D: Risks)
6. True (see Section D: Risks)
7. False (see Section D: Contraindications)
8. False (see Section D: Contraindications)
9. True (see Section D: Precautions)

K. Glossary

This section summarizes important terms used in this information booklet. Please discuss any related questions with your doctor.

Aberration: focusing errors in the eye detectable by wavefront measurements. Examples are nearsightedness and astigmatism (lower-order) and complex errors (higher-order).

Antibiotic Medication: a drug used to treat or prevent infection. Your doctor may prescribe this medication after LASIK surgery.

Anti-inflammatory Medication: a drug that reduces inflammation or the body's reaction to injury or disease. Any eye surgery can cause inflammation. Your doctor may prescribe this medication after LASIK surgery.

Astigmatism: a focusing error that results in blurred distant and/or near vision. The cornea is more curved in some directions than others, and causes light rays to focus at different points inside the eye. Parts of objects appear clearer than other parts.

Autoimmune Disease: a condition in which the body attacks itself and results in inflammation or swelling of parts of the body, such as muscles, joints, and blood vessels. An example is lupus. If you have this type of condition, you should not have LASIK surgery.

Bandage Contact Lens: a soft contact lens placed on the cornea after surgery to cover the area that was treated with the laser.

Cataract: an opacity, or clouding, of the lens inside the eye that can blur vision.

Collagen Vascular Disease: a condition that may result in inflammation or swelling of parts of the body, such as muscles, joints, and blood vessels. An example is rheumatoid arthritis. If you have this type of condition, you should not have LASIK surgery.

Contraindications: any special condition that results in the treatment not being recommended.

Contrast Sensitivity: a measure of the ability of the eye to detect small lightness differences between objects and the background in daylight and in dim light. For example, black lines on a gray background are easier to see than gray lines on a gray background. Objects in daylight are also easier to see than in dim light. Contrast sensitivity testing is a way to determine how well patients can see in poor contrast conditions such as very dim light, rain, snow, and fog.

Conventional LASIK: LASIK surgery that uses an eyeglass prescription to plan the surgery.

Cornea: the clear front layer of the eye. Surgery such as PRK, LASIK and RK reshapes the front surface of the cornea to improve distant vision.

Corneal Flap: a thin slice of tissue on the surface of the cornea made with a microkeratome at the beginning of the LASIK procedure. This flap is folded back before the laser shapes the inner layers of the cornea.

Corneal Swelling: abnormal fluid build-up in the cornea. This condition is usually temporary with no significant effect on vision.

Corneal Wrinkle: temporary appearance of fine white lines in the cornea due to swelling.

CustomCornea® LASIK: LASIK surgery that uses the wavefront to plan the surgery with the LADARVision® 4000 System.

Diopter: a unit of focusing power, used to describe the amount of nearsightedness and astigmatism of an eye. Abbreviated as “D”.

Excimer Laser: a type of laser used in Conventional and CustomCornea® LASIK to remove tissue from the cornea.

Focusing Error: a condition in which your eye forms a blurred image on your retina. Examples are nearsightedness, astigmatism, and higher-order aberrations (complex focusing errors).

Glaucoma: an eye disease usually associated with high eye pressure. Glaucoma damages the optic nerve of the eye and usually causes a progressive loss of vision.

Halos: circular flares or rings of light that may appear around a headlight or other lighted object. This symptom may occur before or after surgery.

Herpes Simplex: a type of viral infection that can recur. This virus typically causes cold sores and/or vesicles to appear on the face or other parts of the body. You should discuss any history of this condition with your doctor before having LASIK surgery.

Herpes Zoster: a type of viral infection that can recur. This condition is a reactivation of the chicken pox virus as an adult. Vesicles appear on only one side of the body. You should discuss any history of this condition with your doctor before having LASIK surgery.

Immunodeficiency Disease: a condition that compromises the body’s ability to heal. An example is acquired immunodeficiency syndrome (AIDS). If you have this type of condition, you should not have LASIK surgery.

Inflammation: the body’s reaction to injury or disease. Eye surgery, such as PRK and LASIK, can cause inflammation.

Keratoconus: a condition of the cornea that results in a change in the shape of the cornea with thinning. If you have this condition, you should not have LASIK surgery.

Laser Assisted In-Situ Keratomileusis (LASIK): a type of eye surgery that uses a microkeratome and a laser to improve vision. The microkeratome creates a thin, hinged flap of tissue on the cornea that is folded back. The laser shapes the tissue under the flap and the flap is put back on the eye so the tissue heals.

Lens: a structure inside the eye that helps to focus light onto the back surface (retina) of the eye.

Low Contrast Visual Acuity: a measure of the sharpness of vision using a 10% low contrast chart with gray letters on a white background. Low contrast acuity testing is another way to determine how well patients can see in poor contrast conditions such as very dim light, rain, snow, and fog.

Microkeratome: a surgical instrument used in LASIK to cut a thin flap of tissue from the front surface of the eye before the laser treatment is applied.

Myopia: a focusing error that results in blurrier vision at distance than near. Myopia is also called nearsightedness.

Nearsightedness: a focusing error that results in blurrier vision at distance than near. Nearsightedness is also called myopia.

Ocular Hypertension: increased eye pressure.

Photorefractive Keratectomy (PRK): a type of eye surgery that uses an excimer laser to reshape the front surface of the eye to improve vision. After the epithelium (outermost layer) of the cornea is first scraped away, the laser removes tissue from the exposed surface. After the surgery, the epithelium grows back.

Radial Keratotomy (RK): a type of eye surgery that changes the shape of the front surface of the eye by making a special pattern of cuts in the cornea to correct nearsightedness.

Retina: the layer of nerve tissue at the back of the eye that captures images, similar to film in a camera, and sends information about those images to the brain. Light must be focused correctly on the retina to form clear images.

Steroid Medication: a drug that reduces inflammation or the body's reaction to injury or disease. Your doctor may prescribe this medication after LASIK surgery for a short time to modify the healing of your eye. If you are taking this medication for a disease condition, you should not have LASIK surgery.

Visual Acuity: a measure of the sharpness of vision using a letter chart.

Wavefront: a measure of the total focusing errors (aberrations) including nearsightedness, astigmatism, and complex focusing errors (higher-order aberrations). Light is projected into your eye and focused on the retina. Part of this light is reflected back out of your eye to form the wavefront.

L. Patient Assistance Information

To be completed by you or your Primary Eye Care Professional as a reference.

Primary Eye Care Professional

Name: _____

Address: _____

Phone: _____

CustomCornea® LASIK Doctor

Name: _____

Address: _____

Phone: _____

Treatment Location

Name: _____

Address: _____

Phone: _____

Laser Manufacturer

| |
|---|
| Alcon, Inc. 2501 Discovery Drive, Suite 500 Orlando, FL 32826 U.S.A. Tel: (877) 523-2784 Fax: (407) 384-1677 |
|---|

M. Index

A

Aberration · 5, *See also* Focusing Error
Definition · 24
Adverse Events · 15
Allergies · 10
Astigmatism · 5
Definition · 6, 24
Autoimmune Disease · 10, 22, 23
Definition · 24

B

Burning · 16, 18

C

Cataract · 9, 20
Definition · 24
Collagen Vascular Disease · 10, 23
Definition · 24
Complications · 15
Contact Lenses
Bandage · 20, 24
Discontinue wear before surgery · 19
Contraindications · 10
Definition · 24
Contrast Sensitivity · 13
Definition · 24
Conventional LASIK · 5
Definition · 24
Cornea · 5
Definition · 25
Flap · 5, 9, 15, 20, 25
Infection · 9, 15
Inflammation · 15
Swelling · 13, 15, 25
Wrinkle · 15, 25
CustomCornea® LASIK · 5
Are You A Good Candidate · 19
Clinical Studies · 8, 13
Definition · 25
How Does It Correct Nearsightedness · 6
Questions To Ask Your Doctor · 21
Self-Test · 22
Summary of Important Information · 23
What Are The Benefits · 8
What Are The Risks · 9
What To Expect During Surgery · 19

D

Diabetes · 10
Diopter · 8
Definition · 25
Double Vision · 16, 18
Double/ghost images · 15
Driving · 8
After Surgery · 19
Dry Eye · 10, 16, 18, 21

E

Excimer Laser · 7
Definition · 25
Eye Tracking · 7, 20
Eyeglass Prescription · 5

F

Feeling of something in the eye · 15, 20
Focusing Error · 5
Definition · 25

G

Glare · 13, 16, 18
Glaucoma · 11, 20
Definition · 25
Gritty Feeling · 16, 18

H

Halos · 16, 18
Definition · 25
Headache · 16, 18
Herpes Simplex · 10
Definition · 25
Herpes Zoster · 10
Definition · 25

I

Immunodeficiency Disease · 10
Definition · 25
Inflammation · 15
Definition · 25
Irregular flap · 9, 15

K

Keratoconus · 10
Definition · 25

L

LADARTracker® System · 7, 20
LADARVision® 4000 System · 5, 7
Clinical Studies · 8, 13
Tracking Eye Movements · 7, 20, 22
Laser Assisted In-Situ Keratomileusis (LASIK) · 5
Definition · 26
Lens · 6
Definition · 26
Light Sensitivity · 13, 16, 18, 20

M

Medication
Amiodarone hydrochloride · 10
Antibiotic · 20, 24
Anti-inflammatory · 13, 20, 24
Isotretinoin · 10
Steroid · 20, 26
Microkeratome · 5, 15, 20
Definition · 26
Myopia · 5, 6, *See also* Nearsightedness
Definition · 26

N

Nearsightedness · 5, 6
Definition · 26
Night driving difficulty · 16, 18

O

Ocular Hypertension · 15, 20
Definition · 26

P

Pain · 13, 15, 16, 18
Photorefractive Keratectomy (PRK) · 5
Definition · 26
Precautions · 11
Pregnancy · 10, 22, 23

R

Radial Keratotomy (RK) · 5
Definition · 26
Reading glasses · 9, 23
Redness · 16, 18
Retina · 6, 7, 15
Definition · 26

T

Tearing · 16, 18

V

Vision
Blurring · 6, 13, 16, 18
Fluctuation · 13, 16, 18
Vision with glasses · 14, 15
Vision without glasses · 8
Visual Acuity · 8, 14
Definition · 26
Low Contrast · 26

W

Warnings · 10
Wavefront · 5
Definition · 26